## Bead One HF Benchtop high frequency fusion instrument







### Bead One HF - Benchtop fusion instrument High frequency induction heating technology


The Bead One HF offers maximum flexibility and accuracy in preparation of glass beads for XRF analysis. Due to its high-frequency induction heating technology the Bead One HF enables complex multi-step fusion processes and shortens preparation time.

### Specially designed for demanding tasks

The Bead One HF offers you extensive functionality and flexibility for demanding analytical tasks. The inductive fusion method allows fast and instantaneous heating-up of the melt to the set temperature. The system operates without inertia and makes it possible to undergo complex preparation cycles requiring extremely precise temperature ramp up. This is useful for preparation procedures of, e.g., ferroalloys depending on several different oxidation and fusion steps of exactly defined temperatures and time lengths.

### Short fusion time and fast ramping up

The Bead One HF reaches temperatures of minimum 1300°C within a short period of time. Accordingly, the Bead One HF is reducing process time and increasing sample throughput. The high fusion and casting dish temperature of the Bead One HF shortens preparation time and improves bead quality for material requiring long-lasting fusion procedures. Moreover, short fusion times lead to low element losses and subsequently improve analysis accuracy and reliability. The crucible is oscillated in a circular movement to improve homogenization of the melt. Additionally, the crucible can be tilt by 45°.



simultaneous heating of crucible and dish



Fusion with seperate pyrometer for crucible and casting dish (optional)



Remarkable features of the Bead One HF are the compact design, easy operation and flexible programmability. The high-frequency induction technology enables complex fusion processes within short period of time.



Easy programm and parameter selection by HMI touch panel



High precision temperature control

# Simultaneous heating of crucible and casting dish

Two separate high frequency generators supply power at the induction coil of the crucible and the casting dish. Since the casting dish is simultaneously heated up the fusion can be casted automatically into the preheated dish. Alternatively, the mold remains in the crucible until solid. The temperature curve of the dish is fully configurable and allows absolutely controlled cooling down subject to analytical requirements. Furthermore, it guarantees that glass buttons are free from reams and seeds. High performance pyrometers for crucible and casting dish ensure absolute temperature constancy. After solidifying, the casting dish is further cooled down by blowing with air.

#### Safety measures

The safety door is locked during the entire fusion and cooling process. Pouring is fully automatically without need for intervention. The Bead One HF has been designed for cold-to-cold operation and the user never comes in contact with hot vessels.

### Easy to operate

Due to its small foot print design the Bead One HF can be set up flexibly at any location. Installation of the machine is very simple, no compressed air is needed. All fusion and cooling parameters can be easily set, modified, stored and monitored using the HMI panel of the machine.

### At a glance

- Benchtop fusion instrument (high frequency induction technology) for preparation of glass disks for XRF analysis
- Highly accurate and precise analysis results for a wide range of different applications
- Optimally suited for complex preparation task due to short ramping up and precise temperature control
- Easy installation and operation, full flexibility in parameter setting, and circular rocking of the crucible for perfectly fused and homogenized beads
- Maximum safety measures including safety door locking and cold-to-cold operation

## **Technical data**

Bead One HF

Colour	Blue/white RAL 5007 / 7035 / 7005
Dimension and Weight LxWxH Weight	800x670x790mm Approx. 140kg
Power supply and cons Voltage Power consumption Protection class	sumption 400 V, 50 Hz, 3/N/PE; CEE 32A 6.5 kVA IP50
Control system and pa Controller Panel Language	nel Siemens S7 Touchpanel 4" German/English
Fusion process Number of definable programms Number of definable steps per programm	<ul><li>16</li><li>Up to 10 parameters:</li><li>Temperature</li><li>Homogenization</li><li>Cooling time</li></ul>
Temperature range crucible Temperature range dish Heating of	385 – 1300 °C 200 – 1200°C
crucible and dish Connections	Simultaneous
Cabel	Terminal
Compressed air Exhaust air stuf Cooling water;	Not required Top of the machine; DN=80mm
inlet/outlet	2x hose nozzle; NW9mm

Cooling water supply Cooling water Inlet pressure Temperature	Min. 2l/min min. 4bar, max. 6bar, cooling water supply directly to the machine 20°C – 30°C	
Bead diameters and cooling		
Diameters	29, 32, 34 and 39mm	
Cooling	Electrical fan	
Input/output		
Crucible, dish, bead	Manual	
Options		
<ul> <li>Temperature measuring of the dish</li> </ul>		
Separate ultrasonic cleaning device		
Separate chiller for water		

ad One HF / 06.2017-E

HERZOG Maschinenfabrik GmbH & Co. KG Auf dem Gehren 1 49086 Osnabrück Germany +49 541 9332-0 Fax +49 541 9332-32

E-Mail info@herzog-maschinenfabrik.de www.herzog-maschinenfabrik.de

HERZOG Automation Corp. 16600 Sprague Road, Suite 400 Cleveland, Ohio 44130 USA

+1 440 891 9777 Fax +1 440 891 9778

E-Mail info@herzogautomation.com www.herzogautomation.com

HERZOG Japan Co., Ltd. 3-7, Komagome 2-chome Toshima-ku Tokio 170-0003, Japan

 +81 3 5907 1771

 Fax
 +81 3 5907 1770

 E-Mail
 nfo@herzog.co.jp

 www.herzog.co.jp

HERZOG (Shanghai) Automation Equipment Co., Ltd. Section A2,2/F, Building 6, No.473, West Fute 1st Road, Waigaoqiao F.T.Z, Shanghai, 200131, P.R. China

+86 21 50375915 Fax +86 21 50375713

E-Mail xc.zeng@herzog-automation.com.cn www.herzog-automation.com.cn

