

TP 20 E

Technical data

Color:	white/blue
Lettering:	English
Instruction manuals:	1 copy, English
Accessories:	1 set of spanners

Dimensions L x W x H:

Machine:	581 x 320 x 537 mm
Machine incl. packaging crate:	800 x 550 x 1,050 mm

Weights:

Machine:	130 kg
Machine incl. packaging crate:	190 kg

Power supply and consumption:

Voltage:	400 V, 50 Hz, 3-phase
Neutral conductor:	Not required
Power consumption:	0.5 kVA
Compressed air supply:	Not required

Electrical switch cabinet:

DDC system:	Siemens LOGO
Type of protection:	IP 54

Process parameters:

Pressure holding time:	To be adjusted individually
Press force:	200 kN
Piston stroke:	32 mm max.

Options:

- Press tool 40 mm diameter for free pressing
- Press tool 40 mm in aluminium cups

- Press tool 40 mm for pressing in steel rings
- Press tool 51.5 mm for pressing in steel rings

The design of the machine complies with the applicable accident prevention and VDE regulations.

We reserve the right to make technical modifications.