# HT 3000 AUTOMATIC CUPWHEEL GRINDER

for Steel and Iron samples







- Compact automatic machine for manual operation and integration into simple automation systems
- possibility to link to HERZOG air tube systems
- little space requirement due to integrated control cabinet
- fast sample handling and good reproducibility of ground sample surface with one grinding wheel
- creation and selection of several grinding programs by the user
- closed, dust-tight and sound-insulated machine housing guaranteeing optimum working security



#### Technical Data

#### Model HT 3000

Colour: blue Lettering: English

Instruction manual: 1-fold

#### **Dimensions**

1050 x 790 x 1600 mm

## Power supply and consumption

Voltage: 400V 50Hz 3Phases Neutral wire: not required Connected load: 10kVA

#### Compressed air supply and consumption

Pressure: min. 5 bar, max. 10 bar Quantity: ca. 750 dm<sup>3</sup>N/Sample

## Water supply and consumption

Pressure: min. 3bar, max. 10 bar Quantity: ca. 2,5 I/Sample

#### **Dedusting connection**

Location: at back

Diameter: OD=120 / ID=115 mm Required dedusting capacity: 15 m³/min. at 2100 Pa Container for dust:

at front

# Control cabinet

PLC-Controler: SIMATIC S 7 Control voltage: 24 V DC

Protection: IP44 Insulation class: B

#### Sample clamping device

Specification: 2 clamping jaws, self-centering and parallel clamping Clamping range - diameter: max. 60 mm Clamping range - thickness: 8...60 mm

#### Grinding parameter

Grinding depth: max. 2 mm Processing time: min. 30 sec,, according to selected program

#### Number of processing programs: 8

#### Sample cooling

Type of cooling: by cooling nozzles

Cooling media: compressed air and water

#### Processable samples

Material: Steel and iron

Shape: without pin, round, oval and dual

thickness sample Hardness: max. 64 HRC Temperature: max: 800° C

#### Sample infeed and outfeed

Infeed method: manual at input position

Location: front at right

Output method: manual or automatic at output

position

#### HERZOG MASCHINENFABRIK GMBH + CO.

Auf dem Gehren 1 D-49086 Osnabrück Tel. 0541 93320 Fax 0541 933232 16600 Sprague Road Cleveland, Ohio 44130, USA Phone 440 891 9777 Fax 440 891 9778

#### HERZOG JAPAN LTD.

5-8, Funado, 3-chome Itabashi-ku, Tokyo 174 Phone 03 5994 0273 Fax 03 3967 2224